





GK3P Antibody

Product Code	CSB-PA614527LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14409
Immunogen	Recombinant Human Glycerol kinase 3 protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerol kinase 3 (GK 3) (Glycerokinase 3) (EC 2.7.1.30) (ATP:glycerol 3-phosphotransferase 3) (Glycerol kinase 3 pseudogene) (Glycerol kinase, testis specific 1), GK3P, GKP3 GKTB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GK3P

Image



Western blot
 All lanes: GK3P antibody at 2µg/ml + Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 61 kDa
 Observed band size: 61 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA614527LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.