



CALR Antibody, Biotin conjugated

Product Code	CSB-PA07867D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27797
Immunogen	Recombinant Human Calreticulin protein (18-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calreticulin (CRP55) (Calregulin) (Endoplasmic reticulum resident protein 60) (ERp60) (HACBP) (grp60), CALR, CRTC
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CALR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.