





ALDH3B2 Antibody

Product Code	CSB-PA001575LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48448
Immunogen	Recombinant Human Aldehyde dehydrogenase family 3 member B2 protein (57-311AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aldehyde dehydrogenase family 3 member B2 (EC 1.2.1.5) (Aldehyde dehydrogenase 8), ALDH3B2, ALDH8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ALDH3B2

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA001575LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA001575LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.