




FGFBP2 Antibody

Product Code	CSB-PA887162LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYJ0
Immunogen	Recombinant Human Fibroblast growth factor-binding protein 2 protein (20-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor-binding protein 2 (FGF-BP2) (FGF-binding protein 2) (FGFBP-2) (37 kDa killer-specific secretory protein) (Ksp37) (HBp17-related protein) (HBp17-RP), FGFBP2, KSP37
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGFBP2

Image



Western blot


All lanes: FGFBP2 antibody at 2µg/ml + Human serum

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA887162LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.