



FGFBP2 Antibody, Biotin conjugated

Product Code	CSB-PA887162LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYJ0
Immunogen	Recombinant Human Fibroblast growth factor-binding protein 2 protein (20-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor-binding protein 2 (FGF-BP2) (FGF-binding protein 2) (FGFBP-2) (37 kDa killer-specific secretory protein) (Ksp37) (HBp17-related protein) (HBp17-RP), FGFBP2, KSP37
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGFBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.