



ATP5F1A Antibody, FITC conjugated

Product Code	CSB-PA002344LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25705
Immunogen	Recombinant Human ATP synthase subunit alpha, mitochondrial protein (44-553AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP synthase subunit alpha, mitochondrial (ATP synthase F1 subunit alpha ATP5A1), ATP5F1A, ATP5A, ATP5A1, ATP5AL2, ATPM
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ATP5F1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.