



MAP2K6 Antibody, Biotin conjugated

Product Code	CSB-PA013415LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52564
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 6 protein (1-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity mitogen-activated protein kinase kinase 6 (MAP kinase kinase 6) (MAPKK 6) (EC 2.7.12.2) (MAPK/ERK kinase 6) (MEK 6) (Stress-activated protein kinase kinase 3) (SAPK kinase 3) (SAPKK-3) (SAPKK3), MAP2K6, MEK6 MKK6 PRKMK6 SKK3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAP2K6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.