



HLA-E Antibody, Biotin conjugated

Product Code	CSB-PA320269LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13747
Immunogen	Recombinant Human HLA class I histocompatibility antigen, alpha chain E protein (33-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class I histocompatibility antigen, alpha chain E (MHC class I antigen E), HLA-E, HLA-6.2 HLAE
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.