




ACBD5 Antibody

Product Code	CSB-PA725893HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T8D3
Immunogen	Recombinant Human Acyl-CoA-binding domain-containing protein 5 protein (1-461AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-CoA-binding domain-containing protein 5, ACBD5, KIAA1996
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ACBD5

Image



Western blot
 All lanes: ACBD5 antibody at 2µg/ml + HeLa whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 61, 55, 59, 47 kDa
 Observed band size: 61 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.