



CCR6 Antibody, FITC conjugated

Product Code	CSB-PA004845EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51684
Immunogen	Recombinant Human C-C chemokine receptor type 6 protein (322-374AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C chemokine receptor type 6 (C-C CKR-6) (CC-CKR-6) (CCR-6) (Chemokine receptor-like 3) (CKR-L3) (DRY6) (G-protein coupled receptor 29) (GPR-CY4) (GPRCY4) (LARC receptor) (CD antigen CD196), CCR6, CKRL3 CMKBR6 GPR29 STRL22
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCR6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.