



G0S2 Antibody, FITC conjugated

Product Code	CSB-PA009114EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27469
Immunogen	Recombinant Human G0/G1 switch protein 2 protein (1-103AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G0/G1 switch protein 2 (G0/G1 switch regulatory protein 2) (Putative lymphocyte G0/G1 switch gene), G0S2
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	G0S2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.