



TNFRSF13C Antibody, Biotin conjugated

Product Code	CSB-PA853495LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RJ3
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 13C protein (1-78AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	B-cell receptor specific for TNFSF13B/TALL1/BAFF/BLyS. Promotes the survival of mature B-cells and the B-cell response.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 13C (B-cell-activating factor receptor) (BAFF receptor) (BAFF-R) (BLyS receptor 3) (CD antigen CD268), TNFRSF13C, BAFFR BR3
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNFRSF13C