



Lhcgr Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA012911LC01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P16235 |
| Immunogen | Recombinant Rat Lutropin-choriogonadotropic hormone receptor protein (27-362AA) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Lutropin-choriogonadotropic hormone receptor (LH/CG-R) (Luteinizing hormone receptor) (LSH-R), Lhcgr |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Neuroscience |
| Target Names | Lhcgr |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |