



Lhcgr Antibody, Biotin conjugated

Product Code	CSB-PA012911LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16235
Immunogen	Recombinant Rat Lutropin-choriogonadotropic hormone receptor protein (27-362AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Relevance	Receptor for lutropin-choriogonadotropic hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lutropin-choriogonadotropic hormone receptor (LH/CG-R) (Luteinizing hormone receptor) (LSH-R), Lhcgr
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Lhcgr