





GPA33 Antibody

Product Code	CSB-PA857459EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99795
Immunogen	Recombinant Human Cell surface A33 antigen protein (22-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface A33 antigen (Glycoprotein A33), GPA33
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	GPA33

Image



Immunofluorescence staining of MCF-7 cells with CSB-PA857459EA01HU at 1:133, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Western Blot

Positive WB detected in: K562 whole cell lysate,
SH-SY5Y whole cell lysate


All lanes: GPA33 antibody at 3.2μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36 kDa

Observed band size: 50 kDa



Immunohistochemistry of paraffin-embedded
human small intestine tissue using CSB-
PA857459EA01HU at dilution of 1:100