



GPA33 Antibody, Biotin conjugated

Product Code	CSB-PA857459ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99795
Immunogen	Recombinant Human Cell surface A33 antigen protein (22-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface A33 antigen (Glycoprotein A33), GPA33
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	GPA33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.