



GPA33 Antibody, Biotin conjugated

Product Code	CSB-PA857459ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99795
Immunogen	Recombinant Human Cell surface A33 antigen protein (22-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface A33 antigen (Glycoprotein A33), GPA33
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	GPA33