



DPH1 Antibody, Biotin conjugated

Product Code	CSB-PA871625LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZG8
Immunogen	Recombinant Human 2-(3-amino-3-carboxypropyl)histidine synthase subunit 1 protein (192-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-(3-amino-3-carboxypropyl)histidine synthase subunit 1 (EC 2.5.1.108) (Diphthamide biosynthesis protein 1) (Diphtheria toxin resistance protein 1) (Ovarian cancer-associated gene 1 protein) (S-adenosyl-L-methionine:L-histidine 3-amino-3-carboxypropyltransferase 1), DPH1, DPH2L DPH2L1 OVCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	DPH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.