




MYO19 Antibody

Product Code	CSB-PA836223LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96H55
Immunogen	Recombinant Human Unconventional myosin-XIX protein (2-770AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unconventional myosin-XIX (Myosin head domain-containing protein 1), MYO19, MYOHD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYO19

Image



Western blot
 All lanes: MYO19 antibody at 2µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 110, 36, 72, 87 kDa
 Observed band size: 110 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.