



DCAF4 Antibody

Product Code	CSB-PA855479LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV16
Immunogen	Recombinant Human DDB1- and CUL4-associated factor 4 protein (1-495AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DDB1- and CUL4-associated factor 4 (WD repeat-containing protein 21A), DCAF4, WDR21 WDR21A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DCAF4

Image



Western blot

All lanes: DCAF4 antibody at 2µg/ml

Lane 1: 293T whole cell lysate

Lane 2: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56, 49, 45, 55 kDa

Observed band size: 56 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.