



MAPK9 Antibody, Biotin conjugated

Product Code	CSB-PA013471LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45984
Immunogen	Recombinant Human Mitogen-activated protein kinase 9 protein (1-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 9 (MAP kinase 9) (MAPK 9) (EC 2.7.11.24) (JNK-55) (Stress-activated protein kinase 1a) (SAPK1a) (Stress-activated protein kinase JNK2) (c-Jun N-terminal kinase 2), MAPK9, JNK2 PRKM9 SAPK1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAPK9