



Tpbg Antibody, HRP conjugated

Product Code	CSB-PA024093LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5PQV5
Immunogen	Recombinant Rat Trophoblast glycoprotein protein (32-361AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trophoblast glycoprotein (5T4 oncofetal trophoblast glycoprotein) (5T4 oncotrophoblast glycoprotein) (Wnt-activated inhibitory factor 1) (WAIF1), Tpbg, 5t4
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Stem Cells
Target Names	Tpbg
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.