



CD86 Antibody

Product Code	CSB-PA004965EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42081
Immunogen	Recombinant Human T-lymphocyte activation antigen CD86 protein (24-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-lymphocyte activation antigen CD86 (Activation B7-2 antigen) (B70) (BU63) (CTLA-4 counter-receptor B7.2) (FUN-1) (CD antigen CD86), CD86, CD28LG2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD86
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.