




Phospho-PRKAB1 (S182) Antibody

Product Code	CSB-PA303153
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y478
Immunogen	Synthesized peptide derived from Human AMPK β 1 around the phosphorylation site of S182.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PRKAB1; AMPK; 5'-AMP-activated protein kinase subunit beta-1; AMPK subunit beta-1; AMPKb
Immunogen Species	Homo sapiens (Human)
Target Names	PRKAB1

Image



Western Blot analysis of COLO205 cells using Phospho-AMPK β 1 (S182) Polyclonal Antibody