





Phospho-INPP5D (Y1021) Antibody

Product Code	CSB-PA903060
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92835
Immunogen	Synthesized peptide derived from Human SHIP-1 around the phosphorylation site of Y1021.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	INPP5D; SHIP; SHIP1; Phosphatidylinositol 3; 4,5-trisphosphate 5-phosphatase 1; Inositol polyphosphate-5-phosphatase of 145 kDa; SIP-145; SH2 domain-containing inositol 5'-phosphatase 1; SH2 domain-containing inositol phosphatase 1; SHIP-1;
Immunogen Species	Homo sapiens (Human)
Target Names	INPP5D

Image



Western Blot analysis of KB cells using Phospho-SHIP-1 (Y1021) Polyclonal Antibody



Western Blot analysis of KB PC-3 cells using Phospho-SHIP-1 (Y1021) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.