





CNP Antibody

Product Code	CSB-PA468752
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09543
Immunogen	Synthesized peptide derived from the Internal region of Human CNPase.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CNP; 2'; 3'-cyclic-nucleotide 3'-phosphodiesterase; CNP; CNPase
Immunogen Species	Homo sapiens (Human)
Target Names	CNP

Image



Western Blot analysis of MCF7 cells using CNPase Polyclonal Antibody



Western Blot analysis of various cells using CNPase Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.