




H1F0 Antibody

Product Code	CSB-PA785142
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07305
Immunogen	Synthesized peptide derived from the Internal region of Human Histone 1.0.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	H1F0; H1FV; Histone H1.0; Histone H1'; Histone H1(0)
Immunogen Species	Homo sapiens (Human)
Target Names	H1F0

Image



Western Blot analysis of A549 cells using Histone 1.0 Polyclonal Antibody