




HNRNPH1/HNRNPH2 Antibody

Product Code	CSB-PA918049
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31943/P55795
Immunogen	Synthesized peptide derived from the Internal region of Human hnRNP H.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HNRNPH1; HNRPH; HNRPH1; Heterogeneous nuclear ribonucleoprotein H; hnRNP H; HNRNPH2; FTP3; HNRPH2; Heterogeneous nuclear ribonucleoprotein H2; hnRNP H2; FTP-3; Heterogeneous nuclear ribonucleoprotein H'; hnRNP H'
Immunogen Species	Homo sapiens (Human)
Target Names	HNRNPH1/HNRNPH2

Image



Western Blot analysis of K562 cells using hnRNP H Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.