




OR51B2 Antibody

Product Code	CSB-PA539310
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5P1
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 51B2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR51B2; OR51B1P; Olfactory receptor 51B2; Odorant receptor HOR5'beta3; Olfactory receptor 51B1
Immunogen Species	Homo sapiens (Human)
Target Names	OR51B2

Image



Western Blot analysis of HT29 cells using Olfactory receptor 51B2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.