




OR51B2 Antibody

Product Code	CSB-PA539310
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5P1
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 51B2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR51B2; OR51B1P; Olfactory receptor 51B2; Odorant receptor HOR5'beta3; Olfactory receptor 51B1
Immunogen Species	Homo sapiens (Human)
Target Names	OR51B2

Image



Western Blot analysis of HT29 cells using Olfactory receptor 51B2 Polyclonal Antibody