




OR52A1 Antibody

Product Code	CSB-PA823392
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKL2
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 52A1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR52A1; Olfactory receptor 52A1; HPFH1OR; Odorant receptor HOR3'beta4; Olfactory receptor OR11-319
Immunogen Species	Homo sapiens (Human)
Target Names	OR52A1

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 52A1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.