




OR52D1 Antibody

Product Code	CSB-PA859541
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H346
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 52D1.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR52D1; Olfactory receptor 52D1; Odorant receptor HOR5'beta14; Olfactory receptor OR11-43
Immunogen Species	Homo sapiens (Human)
Target Names	OR52D1

Image



Western Blot analysis of COS7 cells using Olfactory receptor 52D1 Polyclonal Antibody