




ZC3H7B Antibody

Product Code	CSB-PA000061
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGR2
Immunogen	Synthesized peptide derived from the C-terminal region of Human RoXaN.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZC3H7B; KIAA1031; Zinc finger CCCH domain-containing protein 7B; Rotavirus 'X'-associated non-structural protein; RoXaN
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H7B

Image



Western Blot analysis of 293 cells using RoXaN Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.