




PPP2R5D Antibody

Product Code	CSB-PA245104
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14738
Immunogen	Synthesized peptide derived from the C-terminal region of Human PP2A-B56- δ .
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP2R5D; Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit delta isoform; PP2A B subunit isoform B ¹ -delta; PP2A B subunit isoform B56-delta; PP2A B subunit isoform PR61-delta; PP2A B subunit isoform R5-delta
Immunogen Species	Homo sapiens (Human)
Target Names	PPP2R5D

Image



Western Blot analysis of 293T cells using PP2A-B56- δ Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.