




PRKAA2 Antibody

Product Code	CSB-PA336145
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54646
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	5'-AMP-activated protein kinase catalytic subunit alpha-2
Immunogen Species	Homo sapiens (Human)
Target Names	PRKAA2

Image



Western blot analysis of 1) Hela, 2) 293T, 3) C2C12, 4) 3T3, 5) Rat Heart, 6) Rat Brain using AMPK α 2 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.