





PRKAA1 Antibody

Product Code	CSB-PA801858
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13131
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	5'-AMP-activated protein kinase catalytic subunit alpha-1
Immunogen Species	Homo sapiens (Human)
Target Names	PRKAA1

Image



Immunohistochemical analysis of paraffin-embedded Mouse Skeletal Muscle Tissue using AMPK α 1 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) PC12 using AMPK α 1 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.