




GDAP2 Antibody

Product Code	CSB-PA868350LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXN4
Immunogen	Recombinant Human Ganglioside-induced differentiation-associated protein 2 protein (1-496AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ganglioside-induced differentiation-associated protein 2, GDAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GDAP2

Image



Western blot
All lanes: GDAP2 antibody at 2µg/ml + MCF-7 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 57 kDa
Observed band size: 57 kDa