





KLK5 Antibody

Product Code	CSB-PA897100LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y337
Immunogen	Recombinant Human Kallikrein-5 protein (23-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-5 (EC 3.4.21.-) (Kallikrein-like protein 2) (KLK-L2) (Stratum corneum tryptic enzyme), KLK5, SCTE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KLK5

Image



Western blot
 All lanes: KLK5 antibody at 2µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 32 kDa
 Observed band size: 32 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA897100LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.