



KLK5 Antibody, FITC conjugated

Product Code	CSB-PA897100LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y337
Immunogen	Recombinant Human Kallikrein-5 protein (23-293AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-5 (EC 3.4.21.-) (Kallikrein-like protein 2) (KLK-L2) (Stratum corneum tryptic enzyme), KLK5, SCTE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KLK5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.