




GIMAP7 Antibody

Product Code	CSB-PA818820LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHV1
Immunogen	Recombinant Human GTPase IMAP family member 7 protein (1-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase IMAP family member 7 (Immunity-associated nucleotide 7 protein) (IAN-7), GIMAP7, IAN7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction

Image



Western blot

All lanes: GIMAP7 antibody at 2µg/ml + Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.