




Pklr Antibody

Product Code	CSB-PA018071LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12928
Immunogen	Recombinant Rat Pyruvate kinase PKLR protein (1-574AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate kinase PKLR (EC 2.7.1.40) (L-PK) (Pyruvate kinase isozymes L/R), Pklr
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Pklr

Image



Western blot

All lanes: Pklr antibody at 2µg/ml + Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 63, 59 kDa

Observed band size: 63 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.