



TSPAN31 Antibody, Biotin conjugated

Product Code	CSB-PA622522LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12999
Immunogen	Recombinant Human Tetraspanin-31 protein (94-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-31 (Tspan-31) (Sarcoma-amplified sequence), TSPAN31, SAS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TSPAN31
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.