




TREM1 Antibody

Product Code	CSB-PA024404ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP99
Immunogen	Recombinant Human Triggering receptor expressed on myeloid cells 1 protein (21-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Triggering receptor expressed on myeloid cells 1 (TREM-1) (Triggering receptor expressed on monocytes 1) (CD antigen CD354), TREM1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TREM1

Image



Western blot
 All lanes: TREM1 antibody at 2µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 27, 18, 26 kDa
 Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.