





MMGT1 Antibody

Product Code	CSB-PA014655LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4V1
Immunogen	Recombinant Human Membrane magnesium transporter 1 protein (66-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane magnesium transporter 1 (ER membrane protein complex subunit 5) (Transmembrane protein 32), MMGT1, EMC5 TMEM32
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MMGT1

Image



Immunofluorescence staining of HepG2 cells with CSB-PA014655LA01HU at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

**Western blot**All lanes: MMGT1 antibody at 12 μ g/ml

Lane 1: Hela whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 15, 22 kDa

Observed band size: 15 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.