




CTPS1 Antibody

Product Code	CSB-PA006176LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17812
Immunogen	Recombinant Human CTP synthase 1 protein (300-554AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CTP synthase 1 (EC 6.3.4.2) (CTP synthetase 1) (UTP--ammonia ligase 1), CTPS1, CTPS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CTPS1

Image



Western blot

All lanes: CTPS1 antibody at 6µg/ml

Lane 1: Hela whole cell lysate


Lane 2: 293T whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 67, 42 kDa

Observed band size: 67 kDa



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA006176LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA006176LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.