



CTPS1 Antibody, FITC conjugated

Product Code	CSB-PA006176LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17812
Immunogen	Recombinant Human CTP synthase 1 protein (300-554AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CTP synthase 1 (EC 6.3.4.2) (CTP synthetase 1) (UTP--ammonia ligase 1), CTPS1, CTPS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CTPS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.