



Recombinant Human Interleukin-1 receptor antagonist protein (IL1RN) (Active)

Product Code	CSB-AP001681HU
Abbreviation	Recombinant Human IL1RN protein (Active)
Uniprot No.	P18510
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The ED50 as determined by inhibiting IL-1α- dependent proliferation of murine D10S cells is less than 40 ng/ml, corresponding to a specific activity of >2.5x10 ⁴ IU/mg in the presence of 50 pg/ml rHuIL-1α.
Purity	>96% as determined by SDS-PAGE.
Sequence	RPSGRKSSKM QAFRIWDVNQ KTFYLRNNQL VAGYLQGPV NLEEKIDVVP IEPHALFLGI HGGKMCLSCV KSGDETRLQL EAVNITDLSE NRKQDKRFAF IRSDSGPTTS FESAACPGWF LCTAMEADQP VSLTNMPDEG VMVTKFYFQE DE
Research Area	Immunology
Source	E.coli
Target Names	IL1RN
Expression Region	26-177aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	17.1 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	2139180; 2137201; 1828896; 1385987; 8992991; 1827201; 7629520; 14702039; 15815621; 15489334; 2137200; 2143761; 9514884; 7989776; 7775431; 11991722; 8786086; 16335952; 19494218; 1534997; 8045306; 8175703; 7867645; 9062194; 18987736

Image

**Endotoxin**

Less than 1.0 EU/ μ g as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.