




Recombinant Human Interleukin-31 protein (IL31) (Active)

Product Code	CSB-AP001901HU
Abbreviation	Recombinant Human IL31 protein (Active)
Uniprot No.	Q6EBC2
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, pH7.4
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The specific activity is determined by inducing STAT3 activation using human U-87 MG cells. 5 ng/mL of rHuIL-31 can effectively induce STAT3 activation.
Purity	≥97% as determined by SDS-PAGE.
Sequence	SHTLPVRLLR PSDDVQKIVE ELQSLSKMLL KDVEEEKGVL VSQNYTLPCL SPDAQPPNNI HSPAIRAYLK TIRQLDNKSV IDEIIEHLDK LIFQDAPETN ISVPTDTHC KRFILTISQQ FSECMDLALK SLTSGAQQAT T
Research Area	Immunology
Source	E.coli
Target Names	IL31
Expression Region	24-164aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	15.8 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	15184896; 15489334; 16461143; 17030248

Image





Endotoxin

Less than 1.0 EU/ μ g as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.