




# Recombinant Human Interleukin-36 gamma protein (IL36G) (Active)

<b>Product Code</b>	CSB-AP002031HU
<b>Abbreviation</b>	Recombinant Human IL36G protein (Active)
<b>Uniprot No.</b>	Q9NZH8
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
<b>Product Type</b>	Interleukin
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Biological Activity</b>	Fully biologically active when compared to standard. The specific activity is determined by its binding ability in a functional ELISA. Immobilized rHuIL-36γ at 1 µg/mL can bind recombinant human IL-1 Rrp2 Fc Chimera with a range of 0.15-5 µg/mL.
<b>Purity</b>	≥95% as determined by SDS-PAGE.
<b>Sequence</b>	MRGTPGDADG GGRAVYQSMC KPITGTINDL NQQVWTLQQG NLVAVPRSDS VTPVTVAVIT CKYPEALEQG RGDPIYLGIIQ NPENCLYCEK VGEQPTLQLK EQKIMDLYGQ PEPVKPFLFY RAKTGRTSTL ESVAFPDWF ASSKRDQPII LTSELGKSYN TAFELNIND
<b>Research Area</b>	Immunology
<b>Source</b>	E.coli
<b>Target Names</b>	IL36G
<b>Expression Region</b>	1-169aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	18.7 kDa
<b>Protein Length</b>	Full Length
<b>PubMed ID</b>	10744718; 11466363; 10860666; 11991722; 12975309; 15815621; 15489334; 20870894; 21965679; 23095752; 23147407; 24829417

Image


**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.