




Recombinant Human Interleukin-36 gamma (IL36G) (Active)

Product Code	CSB-AP002041HU
Abbreviation	Recombinant Human IL36G protein (Active)
Uniprot No.	Q9NZH8
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The ED50 as determined by its ability to induce IL-8 secretion by human preadipocytes is less than 10 ng/ml, corresponding to a specific activity of >1x10 ⁵ IU/mg.
Purity	≥95% as determined by SDS-PAGE.
Sequence	SMCKPITGTI NDLNQQVWTL QGQNLVAVPR SDSVTPVTVA VITCKYPEAL EQGRGDPIYL GIQNPEMCLY CEKVGEQPTL QLKEQKIMDL YGQPEPVKPF LFYRAKTGRT STLESVAFPD WFIASSKRQD PIILTSELGK SYNTAFELNI ND
Research Area	Immunology
Source	E.coli
Target Names	IL36G
Expression Region	18-169aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	17.0 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	10744718; 11466363; 10860666; 11991722; 12975309; 15815621; 15489334; 20870894; 21965679; 23095752; 23147407; 24829417

Image


Endotoxin

Less than 1.0 EU/μg as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.